

DATA SHEET

SILICONE PESSARY

MEASURES

FOR.70.51750 diameter 50	FOR.70.51755 diameter 55	FOR.70.51760 diameter 60	FOR.70.51765 diameter 65
FOR.70.51770 diameter 70	FOR.70.51775 diameter 75	FOR.70.51780 diameter 80	FOR.70.51785 diameter 85
FOR.70.51790 diameter 90	FOR.70.51795 diameter 95	FOR.70.517100 diameter 100	/

DESCRIPTION

Device is meant to be used for uterine prolapse.

It is used to substitute the surgical operation on the patient.

This is an invasive device which enters the body across a natural orifice, different from the invasive surgical device, it has not to be connected to an active medical device, and it can be used for a long time (even more than 30 days but not exceed 90 days).

Thickness: 1,2 cm

CLASS OF MEDICAL DEVICE: IIb

GENERAL INFORMATION

Silicone degree/level: GP500

Color: TRANSPARENT

Gives better support and flexibility.

NON-STERILE

INSTRUCTION FOR USE

The Pessary has to be put inside in this way:

1. overcome the vaginal entrance;
2. bend the ellipse ring to help easily the entrance;
3. circle the cervix.

Do not put the device inside in case it falls from the pack, or in case the pack is not intact.


DISPOSAL

Dispose of it at the end of the use, following the local rules about the waste disposal.

No toxic substances are present.

The device is entirely made by silicone parts.

Technical specification:

Producer:	FOR.ME.SA srl 1
Application advice:	Adults  0476
Compliance:	Medical Device Class IIb, comply with Article 12 of Directive 93/42 EEC